



Procedure and Examination Chairs

Top Quality | Versatile | User-friendly | Infinite Positioning



Sonesta manufactures and sells tables and chairs for urology, gynecology, urodynamic diagnostics and upper GI procedures.

Next Generation Sonesta

SONESTA S2



Specifications

Length with extension cushions	193 cm	76"
Length without extension cushions	143 cm	56"
Width	77 cm	30"
Height with seat cushion (min.)	40 cm	16"
Height with seat cushion (max.)	100 cm	39"
Lift capacity	170 kg	375 lbs.
Wheel diameter (locking)	10 cm	4"
Net weight	100 kg	220 lbs.

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All products carry the CE mark

Features:

- Versatility
- Precise positioning
- Two motors offer height and back adjustments
- Easy access for the elderly and handicapped
- Engineered ergonomics
- Infinite adjustments 16" (40 cm) min. to 39" (100 cm) max.
- Lightweight and easy to move
- Supports up to 375 lbs. (170 kg) weight capacity
- UL Certified

Standard Accessories:

2 Arm rests	525-6339
2 Leg rests	525-6336
2 Leg rest holders	525-6338
1 Paper towel holder	525-6326
1 Two motor foot control	524-6353
1 Seat cushion with 170 mm cutout	525-6366-3



We are proud of our Next Generation Sonesta

Our Next Generation Sonesta chairs are the most versatile procedure chairs on the market. The functional design guarantees ultimate comfort and positioning for you and your patient during urologic and gynecologic procedures.

Sonesta is the leading global manufacturer of innovative procedure and examination chairs for urologic and gynecologic diagnostics. Our procedure chairs are manufactured using the highest quality materials to ensure maximum comfort and functionality.



Procedures commonly performed in the Sonesta chairs

- Urodynamics
- Pelvic Muscle Rehabilitation
- Anorectal Manometry
- Pelvic Exam
- Cystoscopy
- Ultrasound
- Vasectomy
- Colposcopy
- Biopsy
- Colonoscopy
- Gynecologic Procedure

Next Generation Sonesta

SONESTA S3



Specifications

Length with extension cushions	193 cm	76"
Length without extension cushions	143 cm	56"
Width	77 cm	30"
Height with seat cushion (min.)	40 cm	16"
Height with seat cushion (max.)	105 cm	41"
Lift capacity	210 kg	465 lbs.
Wheel diameter (Locking)	10 cm	4"
Net weight	125 kg	276 lbs.

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All products carry the CE mark

Features:

- One touch movement locking mechanism for precise positioning
- Versatility
- Three motors offer height, seat and back adjustments
- Easy access for the elderly and handicapped
- Engineered ergonomics
- Infinite adjustments 16" (40 cm) min. to 41" (105 cm) max.
- Lightweight and easy to move
- Supports up to 465 lbs. (210 kg) weight capacity
- UL Certified

Standard Accessories:

- 2 Arm rests 525-6339
- 2 Leg rests 525-6336
- 2 Leg rest holders 525-6338
- 2 Fold up foot rests 525-6331
- 1 Seat cushion without cut out 525-6365
- 1 Seat cushion with 170 mm cutout 525-6366-3
- 1 Paper towel holder 525-6326
- 1 Three motor foot control 524-6357
- 1 Neck pillow 525-6346



Features & benefits of the Next Generation Sonesta chairs

Weight capacity – 375 & 465 lbs.

New soft ergonomic and low profile cushions with improved foam – *To enhance patient comfort and feel.*

Fold up arm rests – *Easy access to the patient from each side: creates a safety precaution when chair is engaged in the table position.*

Height adjustable head rest – *For patient comfort and ease of use.*

Vacuum formed back cover – *Protects and covers the frame and other structural parts, for a calmer and less mechanical feeling. Easier to clean.*

Neck pillow – *Using adjustable magnets for a more premium feel.*

New handlebar – *Supports chair mobility: designed with built-in paper towel holder.*

Footrests – *Newly designed with lever to insure safety for end users. Plastic covered magnet keeps footrest in place when folded up.*

New ergonomic leg rest – *Improved comfort and asymmetric shapes allow leg rest repositioning during procedures if needed.*

Fluoroscopy Procedure Table

SONESTA 6210



Specifications

Length without foot rest/leg extension	156 cm	63"
Length with foot rest/leg extension	196 cm	78"
Width (total)	99 cm	39"
Height with seat cushion (min.)	50 cm	20"
Height with seat cushion (max.)	130 cm	51"
Imaging area without foot rest/leg extension (length)	156 cm	63"
Imaging area with foot rest/leg extension (length)	196 cm	78"
Weight capacity	250 kg	551 lbs.
Weight capacity, standing position		

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All products carry the CE mark

A versatile examination table with C-arm capabilities

The 6210 fluoroscopy procedure table offers the physician standing, seated and supine positioning capabilities. This versatile table easily adjusts from or to standing, seated, and supine positions via a hand held control thus creating a table that is easily accessible to patients of any age. The side-mounted cantilevered design offers the physician anterior/posterior, lateral and oblique views without obstruction. For added convenience, the table is equipped with a battery backup for short-term use in case of power failure or in emergency situations. The superb C-arm access and wide range of positioning possibilities make the 6210 the ideal table for any fluoroscopy procedure.

Features:

- Three motors offer unparalleled positioning and access options.
- Radiolucent carbon fiber seat and back for head-to-toe imaging.
- Four double wheel rolling casters offer unencumbered maneuverability.



Standard Accessories:

1 Hand control	524-6356
1 Plastic funnel w/ drainage	524-6325-9
1 Arm rest holder	524-6338-4
1 Arm rest	524-6021-6
2 Leg rest holders	525-6338
2 Leg rests	525-6336
2 Fold up foot rests	525-6331-10
1 Side rail	524-6329

An ADDvise Group Company

SONESTA

SONESTA MEDICAL AB
 TEGELUDDSVÄGEN 76
 115 28 STOCKHOLM, SWEDEN
 TEL: +46 8 502 572 80
 INFO@SONESTAMEDICAL.SE
 WWW.SONESTAMEDICAL.SE

SONESTA MEDICAL INC.
 2 DEBUSH AVENUE UNIT C3
 MIDDLETON, MA 01949, USA
 TEL: +1 630-519-3450
 INFO@SONESTAMEDICAL.SE
 WWW.SONESTAMEDICAL.SE