

Next Generation Sonesta

SONESTA S2



Features:

- Versatility
- Precise positioning
- Two motors offer height and back adjustments
- Easy access for the elderly and handicapped
- Engineered ergonomics
- Infinite adjustments 16" (40 cm) min. to 39" (100 cm) max.
- Lightweight and easy to move
- Supports up to 375 lbs. (170 kg) weight capacity
- UL Certified

Procedures commonly performed in the Sonesta chairs:

- Urodynamics
- Pelvic Muscle Rehabilitation
- Anorectal Manometry
- Pelvic Exam
- Cystoscopy
- Ultrasound
- Vasectomy
- Colposcopy
- Biopsy
- Colonoscopy
- Gynecologic Procedure



Features & benefits of the Next Generation Sonesta S2 chair

Weight capacity – 170 kg/375 lbs.

New soft ergonomic and low profile cushions with improved foam – *To enhance patient comfort and feel. Easier to clean.*

Fold up arm rests – *Easy access to the patient from each side: creates a safety precaution when chair is engaged in the table position.*

Vacuum formed back cover – *Protects and covers the frame and other structural parts, for a calmer and less mechanical feeling.*

Neck pillow – *Using adjustable magnets for a more premium feel.*

New handlebar – *Supports chair mobility: designed with built-in paper towel holder.*

Footrests – *Newly designed with lever to insure safety for end users. Plastic covered magnet keeps footrest in place when folded up.*

New ergonomic leg rest – *Improved comfort and asymmetric shapes allow leg rest repositioning during procedures if needed.*

Specifications

Length with extension cushions	193 cm	76"
Length without extension cushions	143 cm	56"
Width	77 cm	30"
Height with seat cushion (min.)	40 cm	16"
Height with seat cushion (max.)	100 cm	39"
Lift capacity	170 kg	375 lbs.
Wheel diameter (locking)	10 cm	4"
Net weight	100 kg	220 lbs.

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All products carry the CE mark

Standard Accessories:

2 Arm rests	525-6339
2 Leg rests	525-6336
2 Leg rest holders	525-6338
1 Paper towel holder	525-6326
1 Two motor foot control	524-6353
1 Seat cushion with 170 mm cutout	525-6366-3

