



Specialized in Urology
and Gynecology

SONESTA P2 & P3



SONESTA P2

The Sonesta P2 is an examination chair with two motors for height and back adjustments. It offers the most up-to-date technology, a wide range of accessories and a user-friendly design. The Sonesta P2 creates a pleasant and ergonomic user experience as well as the ultimate patient comfort.

Features & benefits of the Sonesta P2

- Two motors offer height and back adjustments
- Precise positioning and infinite adjustments 40cm (16") to 105cm (41")
- Easily accessible, and easy for the elderly and disabled to enter and exit
- Lightweight and easy to move
- Supports up to 375 lbs. (170 kg) weight capacity



NEWS

- ✓ New electronics with strong motors that run quietly with fast performance
- ✓ Wireless and water proof foot control
- ✓ Three preset and adjustable positionings
- ✓ Electronic trendelenburg
- ✓ 360° rotatable arm rests
- ✓ New mattress color



Specifications, Sonesta P2

Length without extension cushion:	143 cm	76"
Length with extension cushion:	193 cm	56"
Width (total):	69 cm	27"
Height with seat cushion (min.):	40 cm	16"
Height with seat cushion (max.):	105 cm	41"
Patient weight (max):	170 kg	375 lbs
Wheel diameter (2 lockable + 2 directional):	10 cm	4"
Net weight:	80,5 kg	177 lbs

Electric specifications




Power consumption:

Max 3,58A (120V), 1,87A (230V) & 3,9A (100V)

Input voltage: 120 ~ 430VA 50/60 Hz 230V ~ 430VA 50/60 Hz

Output voltage: 24VDC



SONESTA P3

The Sonesta P3 chair is a three motor examination chair that allows for seat, back and lift adjustments. The chair combines sophistication with modern technology, providing more comfort and efficiency during procedures. The Sonesta P3 is the ultimate patient positioning chair with infinite positionings ideal for a variety of urology and gynecology procedures.

Features & benefits of the Sonesta P3

- Three motors offer height, seat and back adjustments
- One touch pelvic locking device
- Precise positioning and infinite adjustments
40 cm (16") to 105 cm (41")
- Easily accessible, and easy for the elderly & disabled to enter and exit
- Lightweight and easy to move
- Supports up to 465 lbs. (210 kg) weight capacity



Specifications, Sonesta P3

Length without extension cushion:	143 cm	76"
Length with extension cushion:	193 cm	56"
Width (total):	69 cm	27"
Height with seat cushion (min.):	40 cm	16"
Height with seat cushion (max.):	105 cm	41"
Patient weight (max):	210 kg	465 lbs
Wheel diameter (2 lockable + 2 directional):	10 cm	4"
Net weight:	83 kg	183 lbs

Electric specifications



Power consumption:

Max 3,58A (120V), 1,87A (230V) & 3,9A (100V)

Input voltage: 120 ~ 430VA 50/60 Hz 230V ~ 430VA 50/60 Hz

Output voltage: 24VDC

NEWS

- ✓ New electronics with strong motors that run quietly with fast performance
- ✓ Wireless and water proof foot control
- ✓ Three preset and adjustable positionings
- ✓ Electronic trendelenburg
- ✓ Automatic pelvic locking device
- ✓ 360° rotatable arm rests
- ✓ New mattress color



Procedures commonly performed in the Sonesta P2 and P3

Urodynamics

- Pelvic Muscle Rehabilitation
- Anorectal Manometry

Cystoscopy

- Pelvic Exam
- Ultrasound

Gynecologic Procedures

- Vasectomy
- Colposcopy
- Biopsy
- Colonoscopy

SONESTA P2

Included accessories

Seat cushion with 170mm cutout	527-120565
Two Motor Hand Control	527-120572

Optional accessories

Two Motor Foot Control	527-120574
Leg Rest Holder (single)	525-6338
Leg Rests (pair)	525-6336
Fold Up Foot Rests (pair)	527-120564
Leg Rest Covers (disposable)	722503
Seat Cushion without cutout	527-120568
Puzzle Piece Extension Cushion	527-120569
Seat Cushion with cutout, 80mm	527-120567
Seat Cushion with cutout, 220mm	527-120566
Seat Covers (disposable)	722504
Neck Pillow	527-120571
Drainage Receptacle	524-6325-1
Funnel Kit	524-6325-5
Funnel & Cysto Pan Kit	524-6325-10
IV Pole	524-6316
Backup Battery	527-120577

SONESTA P3

Included accessories

Seat cushion with 170mm cutout	527-120565
Three Motor Hand Control	527-120573

Optional accessories

Three Motor Foot Control	527-120576
Leg Rest Holder (single)	525-6338
Leg Rests (pair)	525-6336
Fold Up Foot Rests (pair)	527-120564
Leg Rest Covers (disposable)	722503
Seat Cushion without cutout	527-120568
Puzzle Piece Extension Cushion	527-120569
Seat Cushion with cutout, 80mm	527-120567
Seat Cushion with cutout, 220mm	527-120566
Seat Covers (disposable)	722504
Neck Pillow	527-120571
Drainage Receptacle	524-6325-1
Funnel Kit	524-6325-5
Funnel & Cysto Pan Kit	524-6325-10
IV Pole	524-6316
Backup Battery	527-120577



Sonesta Medical AB

Industrivägen 7, 171 48 Solna, Sweden
Tel: +46 8 502 572 80 • info@sonestamedical.se • www.sonestamedical.se

Sonesta Medical Inc.

2 Debush Avenue Unit C3, Middleton, MA 010949, USA
Tel: +1 630-519-3450 • infous@sonestamedical.se • www.sonestamedical.se