

# SONESTA

Specialized in Urology  
and Gynecology

## SONESTA F3



# SONESTA F3

The F3 is a unique procedure table designed for full fluoroscopic imaging of the urinary tract during video Urodynamics procedures, with carbon-fiber support for use with a C-arm. The table has been developed to ensure optimal patient comfort, while focusing on facilitating the positioning of the patient safely and precisely during video fluoroscopy or urology examinations and procedures. The new updated design is optimal for the best hygiene when cleaning.

## NEWS

- ✓ Brand new design
- ✓ Low weight and strong motors with silent and high-speed performance
- ✓ Four pre-set positionings that simplifies examinations
- ✓ Electronic and automatic Trendelenburg
- ✓ Back-up battery
- ✓ Improved safety switch buttons for easier access in emergency situations
- ✓ Hand control with easy touch tactile switches for 12 different positionings of the chair
- ✓ New mattress color
- ✓ Rotatable armrests



## Specifications, Sonesta F3

Length without extension:	156 cm	63"
Length with extension:	194-205 cm	78-81"
Width (total):	104 cm	41"
Height with seat cushion (min.):	51 cm	20"
Height with seat cushion (max.):	138 cm	54"
Patient weight (max):	250 kg	551 lbs
Wheel diameter:	15 cm	5,9"
Net weight:	Approx 252 kg	555 lbs

## Electric specifications



Power consumption:  
Max 5.8A (120V~) 3.0A (230V~)  
Input voltage: 120V ~ 50/60 Hz 230V ~ 50/60 Hz  
Output voltage: 24VDC



# Procedures commonly performed in the Sonesta F3

## Video Urodynamics

- Pelvic Muscle Rehabilitation
- Anorectal Manometry
- Pelvic Exam
- Cystoscopy

## Imaging Examinations

- Ultrasound
- Vasectomy
- Colposcopy

## Fluoroscopy Procedures

- Biopsy
- Colonoscopy
- Gynecologic Procedures

## Features & benefits of the Sonesta F3

- Radiolucent carbon fiber seat and back for head-to-toe imaging
- Three motors enable unparalleled positioning which provides easier access to and visualization of the perineal area
- This versatile chair easily adjusts from/to standing, seated and supine positions, making it easily accessible to all patients
- The versatility allows the patient to remain in one place during the whole urodynamic test and eliminates the risk of catheters falling out or getting misplaced during transfer or ambulation
- The chair can bare patients up to 551 lbs / 250 kg



## SONESTA F3

### Included accessories

Mattress/Cushion with cutout	527-120584
Hand Control	527-120583
Leg Rest Holders (pair)	525-6338
Leg Rests (pair)	525-6336
Arm Rest Left	527-120582
Arm Rest Right	527-120581
Fold Up Foot Rests (pair)	527-120563
Funnel kit (funnel and funnel holder)	524-6325-9
Back-up Battery	527-120577

### Optional accessories

Leg Rest Covers, 150pcs (disposable)	722503
Seat Covers, 20pcs (disposable)	722505
Bolster	722506
Optional Armrest Package	527-120605



**Sonesta Medical AB**

Industrivägen 7, 171 48 Solna, Sweden

Tel: +46 8 502 572 80 • [info@sonestamedical.se](mailto:info@sonestamedical.se) • [www.sonestamedical.se](http://www.sonestamedical.se)

**Sonesta Medical Inc.**

2 Debush Avenue Unit C3, Middleton, MA 010949, USA

Tel: +1 630-519-3450 • [infous@sonestamedical.se](mailto:infous@sonestamedical.se) • [www.sonestamedical.se](http://www.sonestamedical.se)