

SONESTA

Specialized in Urology
and Gynecology

SONESTA P2



SONESTA P2

The Sonesta P2 is an examination chair with two motors for height and back adjustments. It offers the most up-to-date technology, a wide range of accessories and a user-friendly design. The Sonesta P2 creates a pleasant and ergonomic user experience as well as the ultimate patient comfort.


Features & benefits of the Sonesta P2

- Two motors offer height and back adjustments
- Precise positioning and infinite adjustments 40cm (16") to 105cm (41")
- Easily accessible, and easy for the elderly and disabled to enter and exit
- Lightweight and easy to move
- Supports up to 375 lbs. (170 kg) weight capacity



NEWS

- New electronics with strong motors that run quietly with fast performance
- Wireless and water proof foot control
- Three preset and adjustable positionings
- Electronic trendelenburg
- 360° rotatable arm rests
- New mattress color



Specifications, Sonesta P2

Length without extension cushion:	134 cm	53"
Length with extension cushion:	182 cm	72"
Width (total):	69 cm	27"
Height with seat cushion (min.):	40 cm	16"
Height with seat cushion (max.):	105 cm	41"
Patient weight (max):	170 kg	375 lbs
Wheel diameter (2 lockable + 2 directional):	10 cm	4"
Net weight:	80,5 kg	177 lbs

Electric specifications



Power consumption:
Max 3,58A (120V), 1,87A (230V) & 3,9A (100V)
Input voltage: 120 ~ 430VA 50/60 Hz 230V ~ 430VA 50/60 Hz
Output voltage: 24VDC



Procedures commonly performed in the Sonesta P2

[Urodynamics](#)

[Cystoscopy](#)

[Gynecologic Procedures](#)

[Urology](#)

- Pelvic Muscle Rehabilitation
- Anorectal Manometry

- Pelvic Exam
- Ultrasound

- Vasectomy
- Colposcopy

- Biopsy
- Colonoscopy

SONESTA P2

Included accessories

Seat cushion with 170mm cutout	527-120565
Two Motor Hand Control	527-120572

Optional accessories

Two Motor Foot Control	527-120574
Fold Up Foot Rests (pair)	527-120564
Funnel Kit	524-6325-5
Leg Rest Holder	525-6338
Leg Rests (pair)	525-6336
Leg Rest Covers (disposable)	722503
Seat Cushion without cutout	527-120568
Puzzle Piece Extension Cushion	527-120569
Seat Cushion with cutout, 80 mm	527-120567
Seat Cushion with cutout, 220 mm	527-120566
Neck Pillow	527-120571
Drainage Receptacle	524-6325-1
Funnel & Cysto Pan Kit	524-6325-10
IV Pole	524-6316
Backup Battery	527-120577



Sonesta Medical AB

Industrivägen 7, 171 48 Solna, Sweden

Tel: +46 8 502 572 80 • info@sonestamedical.se • www.sonestamedical.se

Sonesta Medical Inc.

2 Debush Avenue Unit C3, Middleton, MA 010949, USA

Tel: +1 630-519-3450 • infous@sonestamedical.se • www.sonestamedical.se