

Product Comparison

SONESTA Chairs

Features	F3	P2	P3
Motorized height adjustment	Yes	Yes	Yes
Motorized seat adjustment	Yes	No	Yes
Motorized reclining back	Yes	Yes	Yes
Seat height (max)	138 cm/54"	105 cm/41"	105 cm/41"
Seat height (min)	51 cm/20"	40 cm/16"	40 cm/16"
Seat tilt range	0°-36°	No	0°-36°
Back tilt range	-15°-90°	-8°-62°	-8°-62°
Weight capacity	250 kg/551 lbs	170 kg/375 lbs	210 kg/465 lbs
Net weight	252 kg/555 lbs	80,5 kg/375 lbs	83 kg/183 lbs
Width	104 cm/41"	69 cm/27"	69 cm/27"
Length without extension cushions	156 cm/63"	134 cm/53"	134 cm/53"
Length with extension cushions	194 cm/78"	182 cm/53"	182 cm/72"
Motor hand control	Yes	Yes	Yes
Motor foot control	No	Option	Option

